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- The BAX® System X5 instrument has a capacity of **32-samples**.
- The BAX® System X5 PCR assays contain **64 tests**
- The computer is a thin 14-inch **laptop**
- The default equipment for the lysis procedure is the compact **DuPont™ Thermal Block**.

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People's  
Republic of  
China AQSIQ

Health  
Canada


Danish  
Veterinary and  
Food Administration

Russian Rospotrebnadzor

Japanese Ministry of Health, Labour  
and Welfare

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
### BAX<sup>®</sup> System Q7 intercalating dye chemistry.

- 1 

PCR reaction tubes contain tablets that are hydrated with prepared test samples. These tablets contain intercalating dye, along with all the reagents, enzymes and primers needed for the PCR reaction. The primers match a unique genetic region of the target's DNA. Their selection determines the specificity of the assay.
- 2 

The PCR process begins with heating the sample solution to denature (separate) double-stranded DNA.
- 3 

A cooling cycle then triggers the highly specific primers to anneal (bind) to a genetic sequence on the single-stranded DNA. As the replicating enzyme begins its copying process, free-floating intercalating dye binds to the newly synthesized double-stranded DNA.
- 4 

Free-floating intercalating dye binds to this region as the DNA re-forms double strands. This chemistry uses a highly specific dye that binds to double-stranded DNA as it accumulates during PCR cycles. When exposed to light, the bound dye emits a measurable fluorescent signal.
- 5 

Successive heating and cooling cycles allow the DNA fragments to replicate (amplify) exponentially within a few hours. The PCR amplification process completes in less than three hours at the end of a cooling cycle with double-stranded DNA. The detection phase begins when samples are exposed to light. This excites the intercalating dye into producing a baseline signal. The temperature is then slowly raised to denature the DNA and measure the resulting drop in signal, which is analyzed as a melt curve.

### BAX<sup>®</sup> System Q7 probe-based chemistry.

- 1 




PCR reaction tubes contain tablets that are hydrated with prepared test samples. These tablets contain highly specific probes, along with all the reagents, enzymes and primers needed for the PCR reaction. The primers and probes match a unique genetic region of the target's DNA. Their selection determines the specificity of the assay.
- 2 

The PCR process begins with heating the sample solution to denature (separate) double-stranded DNA.
- 3 

A cooling cycle then triggers the highly specific primers to anneal (bind) to a genetic sequence on the single-stranded DNA. At the same time, a third primer (probe) also anneals to a region of the sequence. Probes contain a fluorescent reporter dye at one end, whose signal is absorbed by the quencher dye at the other end.
- 4 

As the DNA enzyme extends the primer, it breaks the probe apart. This separates the reporter dye from the quencher, which increases the reporter signal.
- 5 

Successive heating and cooling cycles allow the DNA fragments to replicate (amplify) exponentially. The PCR amplification process completes in about 90 minutes. A separate detection phase is not necessary because signal is measured and analyzed at the end of each cycle during amplification.

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