WATER MICROBIOLOGY TESTING:

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MEDIA PREPARATION

PREPARE M-ENDO SOLUTION

WEIGH DEHYDRATED REAGENT | HYDRATE WITH ETHYL ALCOHOL(95%) | HEAT | ADJUST PH | CAP FLASK | AUTOCLAVE | STORE REFRIGERATED | PREPARE 10/6PDA STOCK SOLUTION | PREPARE MgO₂ STOCK SOLUTION | COMBINE STOCK SOLUTION | CAP FLASK | AUTOCLAVE | RUN CONTROL FOR STERILITY

PREPARE BG18

WEIGH DEHYDRATED REAGENT | HYDRATE WITH H₂O | HEAT | ADJUST PH | ADD DURHAM TUBE | PIPETTE 10ML TO EACH TUBE | CAP TUBES | AUTOCLAVE | RUN POSITIVE CONTROLS | RECORD RESULTS

PREPARE EC MEDIA

WEIGH DEHYDRATED REAGENT | HYDRATE WITH H₂O | HEAT | ADJUST PH | ADD DURHAM TUBE | PIPETTE 10ML TO EACH TUBE | CAP TUBES | AUTOCLAVE | RUN POSITIVE CONTROLS | RECORD RESULTS

RUNNING ROUTINE SAMPLES

LABEL PLATES | PUT PADS ON PLATE | PIPETTE M-ENDO ONTO PLATES | STERILIZE FILTER ASSEMBLY | SET UP FILTER ASSEMBLY | STERILIZE FORCEPS | REMOVE INDIVIDUAL FILTERS | PLACE FILTER ON FILTER ASSEMBLY | RESTERILIZE FORCEPS | SET UP FUNNEL | POUR SAMPLE | TURN ON VACUUM PUMP | RINSE WITH BUFFERED WATER | REMOVE FUNNEL AND STERILIZE | TURN ON UV STERILIZER AND SET TIMER | TRANSFER FILTER TO MENDO PLATES | VERIFY SAMPLE NUMBER | CAP PLATE AND INVERT | INCUBATE 24 HOURS | READ PLATES WITH MICROSCOPE | RECORD RESULTS | RUN 1 BLANK CONTROL EVERY 30 SAMPLES

CONFIRMATION TESTING

COLIFORM CONFIRMATION

SWAB POSITIVE PLATES | LABEL BGLB TUBES | INOCULATE BGLB WITH POSITIVE | CAP TUBES | INCUBATE 24 HOURS | CHECK BGLB FOR GAS AND TURBIDITY | RECORD RESULTS | IF NEGATIVE, INCUBATE 24 HOURS | CHECK BGLB FOR GAS AND TURBIDITY | RECORD RESULTS

FECAL CONFIRMATION

SWAB POSITIVE PLATES | LABEL EC MEDIA TUBES | INOCULATE EC MEDIA WITH POSITIVES | CAP TUBES | INCUBATE 24 HOURS | CHECK FOR GAS AND TURBIDITY | RECORD RESULTS

IDEXX DST®

Presence/Absence

ADD REAGENT TO WATER | INCUBATE VESSEL | CHECK FOR PRESENCE/ABSENCE | RECORD RESULTS

Quantification

ADD REAGENT TO WATER | POUR SAMPLE/PLATE MIXTURE INTO TRAY | SEAL TRAY | INCUBATE TRAY | COUNT YELLOW AND FLUORESCENT WELLS | RECORD RESULTS

There is no comparison.

Call today for a demonstration of how IDEXX can be a smarter solution for all your water testing needs.

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


Worldwide Water Product Approvals/Acceptances

Country	Products Approved ¹	Water Types
Argentina	 ColiIert*	Drinking water
Australia	 ColiIert, ColiIert*-18, Quanti-Tray*, Quanti-Tray*/2000	Drinking water
Belgium	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000, Enterolert*-DW	Drinking water
Bosnia Herzegovina	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Brazil	 ColiIert, ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Canada	 ColiIert, ColiIert-18, Colisure*, Quanti-Tray, Quanti-Tray/2000, SimPlate for HPC	Drinking water
Chile	 ColiIert, ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
China	 ColiIert, Filta-Max*, Quanti-Tray, Quanti-Tray/2000	Drinking water, surface water, ground water, drinking water for poultry and livestock
Colombia	 ColiIert, ColiIert-18, Colisure, Quanti-Tray, Quanti-Tray/2000	Drinking water
Costa Rica	 ColiIert, ColiIert-18, Colisure, Quanti-Tray, Quanti-Tray/2000	Drinking water
Croatia	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water, bottled water
Czech Republic	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Denmark	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
El Salvador	 ColiIert, ColiIert-18, Colisure, Quanti-Tray, Quanti-Tray/2000, SimPlate for HPC	Drinking water
Finland	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water, fresh bathing water
Germany	 ColiIert-18, Enterolert-DW, Quanti-Tray, Quanti-Tray/2000	Drinking water, spas and pools
Greece	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Greenland (Denmark)	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Guatemala	 ColiIert, ColiIert-18, Colisure, Quanti-Tray, Quanti-Tray/2000	Drinking water, source water
Hungary	 ColiIert, ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water, bathing water
Iceland	 ColiIert, ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Ireland	 ColiIert	Drinking water
Italy	 ColiIert-18, Filta-Max, Quanti-Tray, Quanti-Tray/2000	Drinking water, surface water, spas and pools, fresh and marine bathing water

Country	Products Approved ¹	Water Types
Jamaica	 ColiIert, ColiIert-18, Colisure*, Quanti-Tray, Quanti-Tray/2000	Drinking water, source/surface water, wastewater, recreational waters
Japan	 ColiIert	Drinking water
Kuwait	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Latvia	 ColiIert, ColiIert-18, Colisure, Enterolert, Filta-Max, Filta-Max xpress*, Quanti-Tray, Quanti-Tray/2000, SimPlate for HPC, Pseudalert*	Drinking water, bottled water, ambient water, ground water, source/surface water, wastewater, marine water
Lithuania	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Mexico	 ColiIert	Drinking water, bottled water
New Zealand	 ColiIert, ColiIert-18, Colisure, Enterolert*, Quanti-Tray, Quanti-Tray/2000	Drinking water, fresh and marine bathing water
Norway	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water, bottled water
Panama	 ColiIert, ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Philippines	 ColiIert	Drinking water
Poland	 ColiIert-18, Enterolert*-E, Quanti-Disc* YEA, Quanti-Tray, Quanti-Tray/2000	Drinking water, bathing water, source/surface water
Portugal	 ColiIert-18, Quanti-Tray	Drinking water
Slovakia	 ColiIert, ColiIert-18, Colisure, Quanti-Tray, Quanti-Tray/2000	Drinking water, bottled water, ground water, source/surface water, wastewater
South Africa	 ColiIert, Quanti-Tray, Quanti-Tray/2000, Quanti-Disc YEA, SimPlate* for HPC	Drinking water, wastewater effluent
South Korea	 ColiIert	Drinking water
Spain	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000	Drinking water
Sweden	 ColiIert-18, Enterolert-E, Quanti-Tray	Bathing water
Ukraine	 ColiIert-18, Quanti-Tray, Quanti-Tray/2000, Quanti-Disc	Drinking water
United Kingdom	 ColiIert-18, Filta-Max, Quanti-Tray, Quanti-Tray/2000, Quanti-Disc, Enterolert*-DW	Drinking water
United States	 ColiIert, ColiIert-18, Colisure, Enterolert, Filta-Max, Filta-Max xpress*, Quanti-Tray, Quanti-Tray/2000, SimPlate for HPC, Pseudalert*	Drinking water, bottled water, ambient water, ground water, source/surface water, LT2 testing, wastewater, marine water, dairy water (FDA), spas (New Jersey)
Zimbabwe	 ColiIert	Drinking water, ambient water, ground water, wastewater

¡Gracias!

Sonia-O'Donnell@idexx.com

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© IDEXX Laboratories, Inc. 1995